Coapt

# ENDOTINE Midface™-ST 4.5


## Instructions for Use

English	1
Dansk	2
Deutsch	3
Español	5
Français	6
Italiano	8
Nederlands	9
Português	11
Svenska	12
Polski	14
Magyar	15
Ελληνικά	17
Русский	18


CE 0344

1


2


3

4


5


6


7


8

9

10

## ENDOTINE Midface™-ST 4.5

### Available Models

Part No.	Tine Length	Part Name
CFD-060-0197	4.5 mm	ENDOTINE Midface™-ST 4.5

### DESCRIPTION

The ENDOTINE Midface™-ST 4.5 consists of insertion tools and a bioabsorbable implant. The implant consists of a fixation platform attached to an anchoring leash (**Figure 1**). This device and its insertion tools are supplied sterile for single use only.

### INDICATIONS

The ENDOTINE Midface™-ST 4.5 is indicated for use in superperiosteal midface suspension surgery to fixate the cheek subdermis in an elevated position.

### CONTRAINDICATIONS

1. Situations where internal fixation is otherwise contraindicated (e.g., infection).
2. **Patients appearing to have very thin soft tissues of the face in which the implant may be visible or palpable.**
3. Any known allergy or foreign-body sensitivities to plastic biomaterial.

### PACKAGING, LABELING AND STERILIZATION

The ENDOTINE Midface™-ST 4.5 and the supplied insertion tools are pre-sterilized by gamma irradiation and intended for single use only. **Do NOT** resterilize this device or the insertion tools. The device and its insertion tools should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of puncture(s), tampering, water contamination or other damage.

### STORAGE INSTRUCTIONS

1. Store at room temperature (15 to 24°C or 60 to 75°F) in a dry place out of direct sunlight. Do not use beyond the expiration date listed on the outer label.

### INSTRUCTIONS FOR USE

1. The implant device (**Figure 1**) is labeled for orientation including the fixation platform (1A) and anchoring leash (1B).
2. The implant assembly (**Figure 2**) is labeled for orientation including the release mechanism (2A), anchoring leash (2B) and protective cover (2C).
3. The midface dissection from a temporal approach—with or without additional incisions (buccal, periorbital)—is performed using the surgeon's typical approach. In order

to allow the fixation platform to lie over the maxillary antrum after elevation (instead of over the zygoma where it is likely to be visible), the dissection must continue to the inferior maxilla over the buccal sulcus.

4. Insert the implant assembly through the temporal incision into the desired position in the midface (**Figure 3**) with the protective cover oriented towards the skin as illustrated. If significant resistance to the assembly is encountered, do not force it into position. Remove the assembly and extend or enlarge your plane of dissection, with care to protect surrounding nerve structures, in order to provide adequate clearance for proper positioning of the platform at the inferior recess of the dissected tissue.
5. When comfortable with the fixation platform position at the inferior recess of the midface dissection, retract the protective cover by squeezing the release mechanism (**Figure 4**). Apply digital pressure over the cheek to engage the tines in the midface tissue (**Figure 5**), and remove the insertion tools (**Figure 6**).
6. Alternatively, the implant device may be removed from the insertion tools and inserted in a retrograde fashion through an oral incision. (**Figure 7**)
7. Elevate the midfacial tissues to the desired position by applying tension to the anchoring leash (**Figure 8**) exiting the temporal incision. At this point, if the fixation platform appears centered over the zygomatic prominence, undesirable skin irregularities or excessive palpability/sensitivity may result. Beginning the elevation with the device initially more inferior on the midface may avoid this problem.
8. Once the desired position is reached, exercise caution with excessive tissue manipulation or endoscopic visualization, as this may dislodge the fixation platform from the tissue.
9. After elevating the device and midface to the desired position, place suture through the leash (**Figure 9**) and into the temporal fascia where the soft tissue is of adequate thickness.
10. Trim the anchoring leash so that the end lies comfortably within the temporal incision. (**Figure 10**)
11. If adjustment or removal of the implant is desired, one may remove the suture and extract the fixation platform under direct endoscopic guidance. Alternatively, direct visualization of the platform is possible through the buccal incision, if present.
12. The patient should be instructed to exercise caution post-operatively for several weeks to avoid excessive manipulation or contact with the face.
13. Should removal of the device be required in the post-operative period, a small temporal incision to release the leash and an oral approach to remove the platform may offer the most straightforward method.

## PRECAUTIONS

1. Inform the patient of possible adverse effects including palpability of the device prior to resorption, discomfort, visibility of the device under the skin, and surgical risks inherent to the procedure (to include nerve injury, facial musculature dysfunction, infection, scarring, device erosion through the skin, etc.). Some devices may require removal prior to their absorption due to discomfort, infection, reaction or other concerns.
2. Transient local fluid accumulation and/or sinus tract formation may occur. Aspiration may yield traces of implant material and may result in healing of the sinuses without adverse effects.

## WARNINGS

1. DO NOT RESTERILIZE the ENDOTINE Midface™-ST 4.5 or its insertion tools.
2. Do not use if storage temperature exceeds 110°F (43°C) or the temperature indicator on the outer box has turned black.
3. Discard open, unused devices.
4. This device and tools are for SINGLE USE ONLY and MUST NEVER BE RE-USED.
5. Inspect all components for damage before use. Do not attempt to repair damaged components.

## ADVERSE EFFECTS

Material sensitivity/allergic reactions in patients following surgery should be reported. Implantation of foreign materials in tissues can result in histological reactions.

**Caution:** Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

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**Coapt Systems, Inc.**