






Coapt

ENDOTINE TransBleph™

Instructions for Use

English	1
Dansk	2
Deutsch	4
Español	6
Français	7
Italiano	9
Nederlands	11
Português	13
Svenska	14
Polski	16
Magyar	18
Ελληνικά	20
Русский	22

ENDOTINE TransBleph™

Instructions for Use

Available Models

Part Number	Tine Length	Part Name
CFD-080-0177	3.0 mm	ENDOTINE TransBleph™ 3.0
CFD-080-0167	3.5 mm	ENDOTINE TransBleph™ 3.5

DESCRIPTION

The ENDOTINE TransBleph™ consists of insertion tools and a bioabsorbable fixation implant. This device and its insertion tools are for SINGLE USE ONLY and are supplied sterile.

INDICATIONS

The ENDOTINE TransBleph™ is intended for use in subperiosteal browplasty fixation. The ENDOTINE TransBleph™ is specifically indicated for use to fixate the subdermis to the frontal bone.

CONTRAINDICATIONS

1. Patients appearing to have an inadequate frontal bone thickness for the ENDOTINE TransBleph™ bone post which extends to a depth of 2.65 mm. The drilled hole extends to a maximum depth of 3.050 mm. In a series of 54 patients, the inferolateral frontal bone thickness was measured by CT scan. Bone thickness ranged from 4.2 mm to 14.0 mm.
2. Situations where device placement will lead to coverage by eyelid skin instead of brow/forehead tissue.
3. Any plane of dissection other than subperiosteal.
4. Situations where internal fixation is otherwise contraindicated (e.g. infection, sinusitis).
5. Thin, atrophic forehead tissue in which the device might be visible and project excessively.
6. Any known allergy or foreign-body sensitivities to plastic biomaterial.

PACKAGING, LABELING, and STERILIZATION

The ENDOTINE TransBleph™ and its insertion tools are sterilized by gamma irradiation. **Do NOT resterilize this device or the tools. Implants and tools should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of puncture(s), tampering, water contamination or other damage.**

STORAGE INSTRUCTIONS

1. Store at room temperature (15 to 24°C or 60 to 75°F) in a dry place out of direct sunlight. Do not use beyond the expiration date listed on the outer label.

INSTRUCTIONS FOR USE

1. Prior to implant usage, verify intraoperatively that the soft tissue thickness of the forehead is appropriate for the ENDOTINE TransBleph™ implant.
2. The product consists of an implant preloaded on an insertion tool and a dedicated drill bit (**Figure 1**).
3. The supraorbital rim is approached subperiosteally through a standard upper blepharoplasty incision. Care should be exercised to protect the supraorbital and supratrochlear neurovascular bundles. Adequate exposure is required to align the drill bit and insertion tool perpendicular to the bone. If necessary, a small, separate incision in the horizontal lines of expression can be used above the brow for adequate perpendicular access to the frontal bone.
4. The position of the device anchoring hole is planned at the junction of the central third and lateral third of the eyebrow (**Figure 2**). The position of the hole should be an adequate distance above the supraorbital rim to prevent the device from ultimately lying under eyelid skin. Note that when the device is implanted the tines are superior to the post.
5. Use the supplied drill bit to create an anchoring hole perpendicular to the bone surface. The drill bit has reached an adequate depth when the stop collar on the bit is flush with the bone (**Figure 3**). If present, remove any residual loose bone material in the hole after drilling.
6. Remove the ENDOTINE TransBleph™ device from the packaging. Prior to inserting the device into the bone hole, one can pass a double armed suture through the two platform holes in order to later tack down the tissues onto the device. Apply moderate, steady pressure to the insertion tool to seat the implant post in the anchoring hole (**Figure 4**). Rotate the insertion tool counterclockwise to release (**Figure 5**).
7. Elevate the brow tissue to the desired position over the frontal bone.
8. Use digital pressure on the outer surface of the brow to engage brow tissues to the implant tines (**Figure 6**). If the fixation appears inadequate, the suture that was previously passed through the platform holes can be used to fixate subcutaneous tissue or for temporary transcutaneous fixation.
9. Apply gentle pressure to the brow area with surgical tape or elastic fabric to prevent dislodgement in the early post-operative period (**Figure 7**).

PRECAUTIONS

1. The ENDOTINE TransBleph™ device is designed to be seated in the frontal bone. If the frontal sinus or cranial fossa is entered, appropriate consultation should be obtained.
2. Inform the patient of possible adverse effects including palpability and visibility of the device prior to resorption, discomfort, and surgical risks inherent to the procedure (including penetration of the frontal sinus or cranium, nerve injury, facial musculature dysfunction, infection, device erosion through the skin, scarring, inadequate elevation, tissue redescend, etc.). These symptoms and/or wound problems may lead to a decision to remove the device prior to its full resorption.
3. If the device is placed under eyelid skin or thin forehead tissue, it may be visible through the skin and uncomfortable for the patient. Adequate soft tissue coverage is essential.
4. Local cyst and/or sinus formation may occur and may be treated by aspiration or device removal.
5. High speed drilling (approximately 1000 rpm and greater) and repeat drilling may create holes that are wider than intended for a tight grip of the ENDOTINE TransBleph™ insertion post, resulting in inadequate fixation of the post. Single pass, low speed drilling reduces the occurrence of wider holes.
6. The insertion tool is not designed to provide torque and should not be used to twist an implant into position.
7. If a device is removed after anchoring, discard and replace with a new device for optimal performance.

WARNINGS

1. DO NOT RESTERILIZE the ENDOTINE TransBleph™ Device or its insertion tools.
2. Do not cut or reconfigure the ENDOTINE TransBleph™ Device.
3. Do not use if storage temperature exceeds 110°F (43°C) or the temperature indicator on the outer box has turned black.
4. Discard open, unused devices.
5. This device and tools are for SINGLE USE ONLY and MUST NEVER BE REUSED.
6. Inspect all components for damage before use. Do not attempt to repair damaged components or sharpen the drill bit.
7. Incomplete insertion of the device in the frontal bone may result in inadequate anchoring.
8. Contour irregularities of the skin and tenderness may be noticeable overlying the ENDOTINE TransBleph™ device, especially in patients with thinner forehead tissue.

ADVERSE EFFECTS

Material sensitivity/allergic reactions in patients following surgery should be reported. Implantation of foreign materials in tissue can result in histological reactions.

Caution: Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

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Coapt Systems, Inc.